# BRAEBON Q-Snor™ User Guide

For Today's Business of Sleep™

REF 0546

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Guide d'usager du Q-Snor

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- 50 / 60 Hz notch filter: optional
- Sensitivity / Gain:30 uV/mm or 4 000 gain (increase or decrease as needed)

## To use the Q-Snor interface

1. Attach the snore microphone on the naison, cheek or chin using medical grade tape. Alternatively, you may tape the snore microphone on the wall near the sleeping patient. Ensure the snore microphone is pointing towards the nose and mouth.
2. Turn the unit on and connect to the headbox or DC connection. The Q-Snor has one Input (Mic) and 2 Outputs (Volume and Amplified Audio). For the Volume Output use low gain if plugging into a headbox and amplifier; use high gain for a DC connection. The Audio Output can be plugged into the Audio Input of a recording device (PSG acquisition system) to record an amplified Volume Output. Do not connect speakers directly to the Audio Output of the Q-Snor.

When the Q-Snor is turned on, there is a 30-second calibration during which the LED will blink green for 1 second, then red for 1 second, 5 times (10 seconds) with low / high voltage outputs of either 4mV and 10 mV or 0.4V and 1.0V dependent upon the gain switch setting. The LED will then blink red for 10 seconds with a voltage output of 0.4 mV or 0.4 V and green for 10 seconds with a voltage output of either 10 mV or 1.0 V (see Figure 1). After calibration, the LED will blink green once every 6-7 seconds while the Q-Snor is on.

**Warning: Do not connect speakers directly to the audio output of the Q-Snor.**

Table 1. Minimum and maximum volume input and voltage output values.

Gain Switch Setting	Min. Volume & Voltage	Max. Volume & Voltage
Low Gain	40 dB = 4 mV	100 dB = 10 mV
High Gain	40 dB = 0.4 V	100 dB = 1.0 V

Two AA alkaline batteries will last approximately 30 8-hour recordings (i.e., 30 nights). The interface will automatically turn off after ten hours of continuous use. The LED will blink red when it is time to change the batteries.

## Specifications

Operating Temperature: 5°C (41°F)-32°C (90°F)

Storage Temperature: 0°C (32°F)-60°C (140°F)

Operating/Storage Humidity: 10–95% Non-condensing

## Disinfection and Cleaning

**Caution:** Never steam autoclave the devices. Do not immerse the snore microphone in water. Use isopropyl alcohol pads to clean the sensors.

The Q-Snore Microphone is classified as a non-critical item because it only contacts intact skin. You may use isopropyl alcohol pads to gently clean the snore sensor and snore microphone. Do not immerse the snore microphone or Q-Snor interface in any liquid. For additional information, refer to the APIC guidelines for selection and use of disinfectants (American Journal of Infection Control, Vol. 18, number 2, April, 1990, p. 99-117).

## Warranty

BRAEBON MEDICAL CORPORATION warrants to the first consumer that the Q-Snor interface (the "Interface"), when shipped in its original container, will be free from defective workmanship, performance and materials and agrees that it will, at its option, either repair the defect or replace the defective Interface or part thereof at no charge to the purchaser for parts or labor for a time period of twelve months from the date of purchase; the snore microphone (the "Sensor") has a warranty of six months from date of purchase. The warranty described herein shall be the sole and exclusive warranty granted by BRAEBON MEDICAL CORPORATION and shall be the sole and exclusive remedy available to the purchaser. Use of the Interface and / or Sensor constitutes total and complete acceptance of this warranty. Correction of defects, in the manner and for the time period described herein, shall constitute complete fulfillment of all liabilities and responsibilities of BRAEBON MEDICAL CORPORATION to the purchaser with respect to the Interface and / or Sensor and shall constitute full satisfaction of all claims, whether based on contract, negligence, strict liability or otherwise. In no event shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any loss of revenues or damage, direct, incidental, or consequential, including property damage, loss of profit, or personal injury resulting from the use or misuse of, or the inability to use this product. Nor shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any damages or defects in the Interface and / or Sensor which were caused by abuse, misuse, tampering, neglect, incorrect battery type, or repairs or attempted repairs performed by anyone other than an authorized servicer.

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